

	SUBSTRATE COMPATIBILITY	Dok.Nr. 0605-11E
	DIRUNEUTRA	

A Derouging solution composed by a mix of

Diruneutra LIQ
Diruneutra P
Diruclean NS

can be used on the following materials without hesitation:

a) Plastic:

- PP
- PVC
- PE
- PTFE
- PVDF
- EPDM
- Silicone
- PEEK

b) Stainless steel:

- 1.4301
- 1.4401
- 1.4404
- 1.4435
- 1.4571
- 1.4539

c) Glass

Other materials on demand and separate clarification.

Confirmed by	Name	Date	Signature
Ateco Services AG	Marc Vernier	16.12.2009	